

Product Datasheet

Bovine MUC1 ELISA Kit orb403317

Description	Bovine Mucin-1(MUC1) ELISA kit
Reactivity	Bovine
Range	0.312 U/ml-20 U/ml
Sensitivity	0.078 U/ml.
Tested Applications	ELISA
Target	MUC1
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm