

## Product Datasheet

### Bovine TGFB2 ELISA Kit orb403325

<b>Description</b>	Bovine transforming growth factors $\beta$ 2, TGF $\beta$ 2 ELISA Kit
<b>Reactivity</b>	Bovine
<b>Range</b>	6.25 ng/ml-100 ng/ml
<b>Sensitivity</b>	6.25 ng/ml
<b>Tested Applications</b>	ELISA
<b>Target</b>	TGFB2
<b>Note</b>	For research use only.
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	<ul style="list-style-type: none"><li>• <b>Sample Type:</b> serum, plasma, tissue homogenates</li><li>• <b>Assay Time:</b> 1-5 hours</li><li>• <b>Sample Volume:</b> 50-100 <math>\mu</math>l</li><li>• <b>Detection Wavelength:</b> 450 nm</li></ul>