

Product Datasheet

Canine CALP ELISA Kit orb403328

Description	Canine calprotectin (CALP) ELISA kit
Reactivity	Canine
Range	6.25 ng/ml-400 ng/ml
Sensitivity	1.56 ng/ml
Tested Applications	ELISA
Target	CALP
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates, cell lysates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm