

## Product Datasheet

### Chicken CD4 ELISA Kit orb403353

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Chicken Cluster of Differentiation 4(CD4) ELISA Kit   |
| <b>Reactivity</b>          | Gallus  |
| <b>Range</b>               | 3.125 ng/ml-800 ng/ml   |
| <b>Sensitivity</b>         | 3.125ng/ml  |
| <b>Tested Applications</b> | ELISA   |
| <b>Target</b>              | CD4   |
| <b>Note</b>                | For research use only.  |
| <b>Expiration Date</b>     | Please enquire.   |
| <b>Application Notes</b>   | <ul style="list-style-type: none"><li>• <b>Sample Type:</b> serum, plasma, tissue homogenates</li><li>• <b>Assay Time:</b> 1-5 hours</li><li>• <b>Sample Volume:</b> 50-100 <math>\mu</math>l</li><li>• <b>Detection Wavelength:</b> 450 nm</li></ul> |