

Product Datasheet

Chicken CD8 ELISA Kit (orb403354)

Description	Chicken Cluster of differentiation 8,CD8ELISA Kit
Reactivity	Gallus
Range	3.125 U/mL-800 U/mL
Sensitivity	3.125 U/mL
Tested Applications	ELISA
Target	CD8
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q6QR64
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm