
Product Datasheet

Chicken IFN gamma ELISA Kit (orb403361)

Description	Chicken Interferon gamma; ,IFN-gamma; ELISA Kit
Reactivity	Gallus
Range	3.125 pg/mL-800 pg/mL
Sensitivity	3.125 pg/mL
Tested Applications	ELISA
Target	IFNG
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P49708
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm