

Product Datasheet

Chicken IL-12/70 ELISA Kit orb403364

Description	Chicken Interleukin 12,IL-12/70 ELISA Kit
Reactivity	Gallus
Range	0.78 pg/ml-200 pg/ml
Sensitivity	0.78 pg/ml
Tested Applications	ELISA
Target	IL-12/70
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma and tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm