

Product Datasheet

Chicken IL-17 ELISA Kit orb403366

Description	Chicken Interleukin 17,IL-17 ELISA Kit
Reactivity	Gallus
Range	150 pg/ml-6000 pg/ml
Sensitivity	150 pg/ml
Tested Applications	ELISA
Target	IL-17
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm