

## Product Datasheet

### Chicken GPX-8 ELISA Kit orb403386

<b>Description</b>	Chicken Glutathione Peroxidase 8 (GPX-8) ELISA Kit
<b>Reactivity</b>	Gallus
<b>Range</b>	31.25 mU/mL-2000 mU/mL
<b>Sensitivity</b>	7.8 mU/mL
<b>Tested Applications</b>	ELISA
<b>Target</b>	GPX-8
<b>Note</b>	For research use only.
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<b>R4GGB2</b>
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm