

Product Datasheet

Chicken streptolysin O(ASO) antibody (IgG) ELISA Kit orb403396

Description	Chicken anti-streptolysin O(ASO) antibody (IgG) ELISA Kit
Reactivity	Gallus
Range	Qualitative
Tested Applications	ELISA
Target	streptolysin O(ASO) antibody (IgG)
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm