

Product Datasheet

Chicken growth hormone ELISA Kit orb403416

Description	Chicken growth hormone,GH ELISA KIT
Reactivity	Gallus
Range	625 pg/mL-10000 pg/mL
Sensitivity	625 pg/mL
Tested Applications	ELISA
Target	GH
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P08998
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm