

Product Datasheet

Chicken growth hormone ELISA Kit orb403416

Description	Chicken growth hormone,GH ELISA KIT
Reactivity	Gallus
Range	625 pg/ml-10000 pg/ml
Sensitivity	625 pg/ml
Tested Applications	ELISA
Target	growth hormone
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm