

Product Datasheet

Chicken INS ELISA Kit orb403425

Description	Chicken Insulin ELISA kit
Reactivity	Gallus
Range	3.12 µIU/ml-50 µIU/ml
Sensitivity	0.59 µIU/ml
Tested Applications	ELISA
Target	INS
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm