

Product Datasheet

Chicken TGF beta 1 ELISA Kit orb403434

Description	Chicken Transforming Growth factor β 1, TGF- β 1 ELISA kit
Reactivity	Gallus
Range	0.312 ng/ml-20 ng/ml
Sensitivity	0.078ng/ml
Tested Applications	ELISA
Target	TGF beta 1
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm