

Product Datasheet

Animal IGF-1 ELISA Kit (orb403437)

Description	Deer Insulin-like growth factors 1,IGF-1 ELISA Kit
Reactivity	Animal
Range	1 ng/mL-200 ng/mL
Sensitivity	0.5 ng/mL
Tested Applications	ELISA
Target	IGF-1
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	H9BDM1
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm