

Product Datasheet

Dog IL2RA ELISA Kit orb403471

Description	Dog Interleukin-2 receptor subunit alpha(IL2RA) ELISA kit
Reactivity	Canine
Range	62.5 pg/ml-4000 pg/ml
Sensitivity	15.6 pg/ml.
Tested Applications	ELISA
Target	IL2RA
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm