

Product Datasheet

Dog IL33 ELISA Kit orb403474

Description	Dog Interleukin-33(IL33) ELISA kit
Reactivity	Canine
Range	23.44 pg/ml-1500 pg/ml
Sensitivity	5.86 pg/ml
Tested Applications	ELISA
Target	IL33
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm