

Product Datasheet

Dog HCRT ELISA Kit orb403483

Description	Dog Orexin(HCRT) ELISA kit
Reactivity	Canine
Range	1.87 pg/ml-120 pg/ml
Sensitivity	0.47 pg/ml
Tested Applications	ELISA
Target	HCRT
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm