
Product Datasheet

Dog HCRT ELISA Kit (orb403483)

Description	Dog Orexin(HCRT) ELISA kit
Reactivity	Canine
Range	1.87 pg/mL-120 pg/mL
Sensitivity	0.47 pg/mL
Tested Applications	ELISA
Target	HCRT
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9GLF6
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm