

Product Datasheet

Dog NT-proBNP ELISA Kit orb403505

Description	Canine N-terminal pro-brain natriuretic peptide(NT-proBNP)ELISA Kit
Reactivity	Canine
Range	0.312 ng/ml-20 ng/ml
Sensitivity	0.078 ng/ml
Tested Applications	ELISA
Target	NT-proBNP
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm