

Product Datasheet

Dog S100B ELISA Kit orb403507

Description	Dog S100 calcium binding protein B (S100B) ELISA kit
Reactivity	Canine
Range	39 pg/mL-2500 pg/mL
Sensitivity	9.8 pg/mL
Tested Applications	ELISA
Target	S100B
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates_x000D_
Uniprot ID	E2QTP3
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm