

---

## Product Datasheet

### Dog IFN- $\alpha$ 1/2 ELISA Kit (orb403509)

<b>Description</b>	Canine Interferon alpha-1/2 ELISA kit
<b>Reactivity</b>	Canine
<b>Range</b>	12.5 pg/mL-800 pg/mL
<b>Sensitivity</b>	3.12 pg/mL
<b>Tested Applications</b>	ELISA
<b>Target</b>	IFN-A1/2
<b>Note</b>	For research use only.
<b>Sample Types</b>	serum, plasma, cell culture supernates, tissue homogenates
<b>Uniprot ID</b>	<b>P81255</b>
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	• Sample Type: serum, plasma, cell culture supernates, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 $\mu$ l • Detection Wavelength: 450 nm