

Product Datasheet

Dog IFN- α 1/2 ELISA Kit orb403509

Description	Canine Interferon alpha-1/2 ELISA kit
Reactivity	Canine
Range	12.5 pg/ml-800 pg/ml
Sensitivity	3.12 pg/ml
Tested Applications	ELISA
Target	IFN- α 1/2
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, cell culture supernates, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm