

Product Datasheet

Dog IL17 ELISA Kit orb403517

Description	Canine Interleukin 17(IL-17)ELISA Kit
Reactivity	Canine
Range	31.25 pg/ml-2000 pg/ml
Sensitivity	7.81 pg/ml.
Tested Applications	ELISA
Target	IL17
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm