

Product Datasheet

Dog C5a ELISA Kit [Out of Stock] orb403529

Description	Dog complement fragment 5a (C5a) ELISA kit
Reactivity	Canine
Range	7.8 ng/ml-500 ng/ml
Sensitivity	1.95 ng/ml
Tested Applications	ELISA
Target	C5a
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-3 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450nm