

Product Datasheet

Dog SOD2 ELISA Kit [Out of Stock] orb403531

Description	Dog Superoxide dismutase [Mn], mitochondrial(SOD2) ELISA kit
Reactivity	Canine
Range	1.25 ng/ml-80 ng/ml
Sensitivity	0.312 ng/ml.
Tested Applications	ELISA
Target	SOD2
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm