

Product Datasheet

Dog Luteinizing Hormone ELISA Kit orb403549

Description	Canine Luteinizing Hormone,LH ELISA Kit
Reactivity	Canine
Range	0.75 mIU/mL-60 mIU/mL
Sensitivity	0.94 mIU/mL
Tested Applications	ELISA
Target	LH
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P18842
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma and tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm