

Product Datasheet

Dog MMP3 ELISA Kit (orb403555)

Description	Dog Matrix metalloproteinase 3, MMP-3 ELISA kit
Reactivity	Canine
Range	1.56 ng/mL-100 ng/mL
Sensitivity	0.39 ng/mL
Tested Applications	ELISA
Target	MMP3
Note	For research use only.
Sample Types	serum, plasma, cell culture supernates, urine, tissue homogenates
Uniprot ID	Q6Y4Q5
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, cell culture supernates, urine, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm