

Product Datasheet

Dog PTH ELISA Kit orb403557

Description	Canine Parathyroid hormone(PTH) ELISA Kit
Reactivity	Canine
Range	15.6 pg/mL-1000 pg/mL
Sensitivity	3.9 pg/mL
Tested Applications	ELISA
Target	PTH
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P52212
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm