

Product Datasheet

Dog IgA ELISA Kit orb403568

| | |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description | Canine Immunoglobulin A,IgA ELISA Kit |
| Reactivity | Canine |
| Tested Applications | ELISA |
| Target | IgA |
| Note | For research use only. |
| Expiration Date | Please enquire. |
| Application Notes | <ul style="list-style-type: none">• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm |