

Product Datasheet

Dog MPO_ Myeloperoxidase ELISA Kit orb403589

Description	Dog Myeloperoxidase (MPO) ELISA kit
Reactivity	Canine
Range	31.25 ng/ml-2000 ng/ml
Sensitivity	7.81 ng/ml
Tested Applications	ELISA
Target	MPO_ Myeloperoxidase
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, cell culture supernates and tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm