

Product Datasheet

Dog TARC/CCL17 ELISA Kit orb403599

Description	Canine thymus activation regulated chemokine(TARC/CCL17)ELISA Kit
Reactivity	Canine
Range	23.5 pg/mL-1500 pg/mL
Sensitivity	5.86 pg/mL
Tested Applications	ELISA
Target	TARC/CCL17
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q95N01
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm