

Product Datasheet

Dog ANG-II ELISA Kit orb403602

Description	Dog angiotensin II (ANG-II) ELISA Kit
Reactivity	Canine
Range	6.25 pg/mL-400 pg/mL
Sensitivity	1.56 pg/mL
Tested Applications	ELISA
Target	ANG-II
Note	For research use only.
Sample Types	serum, plasma
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm