

Product Datasheet

Dog VWF ELISA Kit orb403603

Description	Dog von Willebrand Factor,vWF ELISA Kit
Reactivity	Canine
Range	0.625 ng/ml-40 ng/ml
Sensitivity	0.32 ng/ml
Tested Applications	ELISA
Target	VWF
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma and tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm