

Product Datasheet

Dog VWF ELISA Kit orb403603

Description	Dog von Willebrand Factor,vWF ELISA Kit
Reactivity	Canine
Range	0.625 ng/mL-40 ng/mL
Sensitivity	0.32 ng/mL
Tested Applications	ELISA
Target	VWF
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q28295
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma and tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm