

Product Datasheet

Dog PNLIP ELISA Kit orb403615

Description	Dog Pancreatic triacylglycerol lipase(PNLIP) ELISA kit
Reactivity	Canine
Range	0.45 U/ml-30 U/ml
Sensitivity	0.11 U/ml
Tested Applications	ELISA
Target	PNLIP
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm