

## Product Datasheet

### Dog Complement 4 ELISA Kit orb403624

<b>Description</b>	Dog Complement 4(C4) ELISA kit
<b>Reactivity</b>	Canine
<b>Range</b>	31.25 ng/mL-2000 ng/mL
<b>Sensitivity</b>	7.8 ng/mL
<b>Tested Applications</b>	ELISA
<b>Target</b>	C4
<b>Note</b>	For research use only.
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-3 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450nm