

Product Datasheet

Dog CIC ELISA Kit orb403626

Description	Dog circulating immune complex (CIC) ELISA kit
Reactivity	Canine
Range	2.34 ng/ml-1200 ng/ml
Sensitivity	2.34 ng/ml
Tested Applications	ELISA
Target	CIC
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-3 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450nm