

Product Datasheet

Bird IFN- γ ELISA Kit orb403627

Description	Duck Interferon gamma; ,IFN-gamma; ELISA Kit
Reactivity	Aves
Range	15.6 pg/mL-1000 pg/mL
Sensitivity	3.9 pg/mL
Tested Applications	ELISA
Target	IFN-
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9YGB9
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 μ l • Detection Wavelength: 450 nm