

Product Datasheet

Fish IL- 1 β ELISA Kit orb403635

Description	Fish Interleukin 1 β (IL- 1 β) ELISA Kit
Reactivity	Fish
Range	3.12 pg/ml -200 pg/ml
Sensitivity	0.78 pg/ml
Tested Applications	ELISA
Target	IL- 1 β
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm