

## Product Datasheet

### Fish IL8/CXCL8 ELISA Kit orb403637

<b>Description</b>	Fish Interleukin-8,IL-8 ELISA Kit
<b>Reactivity</b>	Fish
<b>Range</b>	31.25 pg/mL-500 pg/mL
<b>Sensitivity</b>	31.25 pg/mL
<b>Tested Applications</b>	ELISA
<b>Target</b>	IL8/CXCL8
<b>Note</b>	For research use only.
<b>Sample Types</b>	serum, plasma
<b>Uniprot ID</b>	<b>Q6F4N2</b>
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	• Sample Type: serum, plasma • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm