

Product Datasheet

Fish ACTH ELISA Kit orb403640

Description	Fish adrenocorticotrop hormone (ACTH) ELISA kit
Reactivity	Fish
Range	75 pg/ml-1200 pg/ml
Sensitivity	75 pg/ml
Tested Applications	ELISA
Target	ACTH
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm