

## Product Datasheet

### Fish MT ELISA Kit orb403647

<b>Description</b>	Fish metallothionein(MT)ELISA kit
<b>Reactivity</b>	Fish
<b>Range</b>	25 pg/ml-400 pg/ml
<b>Sensitivity</b>	12.5 pg/ml
<b>Tested Applications</b>	ELISA
<b>Target</b>	MT
<b>Note</b>	For research use only.
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	<ul style="list-style-type: none"><li>• <b>Sample Type:</b> serum, plasma</li><li>• <b>Assay Time:</b> 1-5 hours</li></ul>