

Product Datasheet

Fish TNF- α ELISA Kit orb403662

Description	Fish Tumor Necrosis Factor α (TNF- α) ELISA KIT
Reactivity	Fish
Range	125 pg/ml-5000 pg/ml
Sensitivity	125 pg/ml
Tested Applications	ELISA
Target	TNF- α
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm