

Product Datasheet

Fish TNF- α ELISA Kit (orb403662)

Description	Fish Tumor Necrosis Factor α (TNF- α) ELISA KIT
Reactivity	Fish
Range	125 pg/mL-5000 pg/mL
Sensitivity	125 pg/mL
Tested Applications	ELISA
Target	TNF- α
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q8JFG3
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 μ l • Detection Wavelength: 450 nm