

Product Datasheet

Goat CSN2 ELISA Kit orb403668

Description	Goat Beta-casein(CSN2) ELISA kit
Reactivity	Goat
Range	0.007 mg/ml -15 mg/ml
Sensitivity	0.007 mg/ml
Tested Applications	ELISA
Target	CSN2
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: breast milk• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm