

Product Datasheet

Goat Myoglobin ELISA Kit [Out of Stock] orb403673

Description	Goat Myoglobin(MB) ELISA kit
Reactivity	Goat
Range	31.25 pg/ml-2000 pg/ml
Sensitivity	7.82 pg/ml
Tested Applications	ELISA
Target	Myoglobin
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates.• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm