

Product Datasheet

Goat 25-Dihydroxyvitamin ELISA Kit orb403677

Description	Goat 1,25-Dihydroxyvitamin D3(1,25 DHVD3) ELISA kit
Reactivity	Goat
Tested Applications	ELISA
Target	25-Dihydroxyvitamin
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm