

Product Datasheet

Goat IL12B ELISA Kit orb403684

Description	Goat interleukin 12,IL-12/P40 ELISA KIT
Reactivity	Goat
Range	9.38 pg/ml-600 pg/ml
Sensitivity	2.34pg/ml
Tested Applications	ELISA
Target	IL12B
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates, cell lysates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm