

Product Datasheet

Goat IL2 ELISA Kit orb403688

Description	Goat Interleukin 2(IL-2) ELISA Kit
Reactivity	Goat
Range	15.6 pg/ml-1000 pg/ml
Sensitivity	3.9 pg/ml
Tested Applications	ELISA
Target	IL2
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm